


# Accessories



**1002-1008**



**Test Tube Rack**  
made of stainless steel,  
with 20 openings of 18 mm Ø,  
for test tubes of 16/17 mm Ø,  
max. height 185 mm

**Order No. 1920**

Capacities:

Water Bath 1002	4 racks =	80 test tubes
Water Bath 1003 / 1008	8 racks =	160 test tubes
Water Bath 1004	12 racks =	240 test tubes
Water Bath 1005	9 racks =	180 test tubes




**Test Tube Rack**  
made of stainless steel,  
with 5 openings of 31 mm Ø,  
max. height 185 mm

**Order No. 1921**

Capacities:

Water Bath 1002	4 racks =	20 test tubes
Water Bath 1003 / 1008	8 racks =	40 test tubes
Water Bath 1004	12 racks =	60 test tubes
Water Bath 1005	9 racks =	45 test tubes



**Test Tube Rack**  
made of stainless steel,  
with 20 openings of 13 mm Ø,  
for test tubes of 12 mm Ø,  
max. height 185 mm

**Order No. 1922**

Capacities:

Water Bath 1002	4 racks =	80 test tubes
Water Bath 1003 / 1008	8 racks =	160 test tubes
Water Bath 1004	12 racks =	240 test tubes
Water Bath 1005	9 racks =	180 test tubes



**Milk Bottle Rack**  
made of stainless steel,  
with 12 openings of 56 mm Ø,  
for baby milk bottles

**Order No. 1942**

Capacities:

Water Bath 1002	1 rack =	12 bottles
Water Bath 1003 / 1008	2 racks =	24 bottles
Water Bath 1004	3 racks =	36 bottles
Water Bath 1005	2 racks =	24 bottles




**Adjustable  
Water Level Regulator**  
to keep the water level constant  
and to cool Water Baths 1002-  
1008 (approx. 3 °C above tap  
water temperature)

**Order No. 1919**



**1005**



**Rack for Water Bath 1005**  
made of stainless steel,  
especially suitable for warming  
medical hot packs.

**Order No. 1923**

Capacities:

Water Bath 1005	up to 8 hot packs
-----------------	-------------------