


# 1023

## Indispensable


For gentle steaming work with Erlenmeyer flasks, glasses, etc., GFL Water Baths models Steam Bath 1023 and Water Baths for fume hoods 1031 and 1032 have become indispensable in daily practice.

### Accessories



**Test Tube Rack**  
made of stainless steel,  
with 100 openings of  $\varnothing$  18 mm,  
for test tubes

Order No. 1933



**Support Rod**  
made of stainless steel,  
length 316 mm,  $\varnothing$  12 mm

Order No. 1985

GFL Water Baths convince through their varied ranges of application. For instance, the Steam Bath's square cover frame (W 265 mm x L 265 mm) is removable. Its 9-part set of rings, made of heat-resistant plastic material can be parted. Its diameter can be varied in approx. 20 mm steps (min. 32.5 / max. 173.5 mm).

1023 Steam Bath, 7 litres



### Technical data

Model Order No.	Interior dimensions			Exterior dimensions			Capacity approx. litres	Electrical connection*	Weight approx.		Packing volume approx m <sup>3</sup>
	Width mm	Depth mm	Working height <sup>♦</sup> mm	Width mm	Depth mm	Height mm			net kg	gross kg	
1023	240	240	120	342	400	180	7	230 V / 50...60 Hz / 1.0 kW	9	12	0.10
	♦ perforated floor to upper rim			* Other voltages available on request							