




Heating mantle | beaker

for beaker, integrated magnetic stirrer and controller, up to 450°C



1.000 ml beaker heating mantle

WHM for beaker, 100°C - 450°C with built-in magnetic stirrer integrated controller

Ideal for:

- for beaker 100 - 5.000 ml

Features:

- built-in stirring and temperature controller
- 50 ~ 750 rpm
- full choice of vessel, capacity, housing and control are offered
- nickel chrome heating element
- coil type: acid-resistant, long life time
- all kinds of beaker can be used (Duran®/Kimax/Pyrex/Iwaki)
- CE certified and unique serial number for tracing

Controller:

- built-in *stirring & temperature controller*



2x stirrer bars included

Application (beaker brand)	Beaker capacity (ml)	Beaker Ø (mm)	Watts (W)	Inside depth (mm)	Mantle outside Ø x H mm	Support rod Ø x H mm	Order number
							with 230V
common-use	100	50 - 52	80	45	Ø 145 x 140	Ø 13 x 440	DH.WHM121601
common-use	250	68 - 70	140	70	Ø 190 x 165	Ø 13 x 440	DH.WHM121602
for DU.	400	79 - 81	240	70	Ø 190 x 165	Ø 13 x 440	DH.WHM121603
for Iw./Ki./Py.	400	76 - 78	240	70	Ø 190 x 165	Ø 13 x 440	DH.WHM121604
common-use	600	88 - 90	240	72	Ø 190 x 165	Ø 13 x 440	DH.WHM121605
common-use	1,000	105 - 108	350	90	Ø 190 x 200	Ø 13 x 440	DH.WHM121606
common-use	2,000	130 - 132	450	140	Ø 270 x 235	Ø 13 x 660	DH.WHM121607
for DU./Iw.	3,000	151 - 153	550	140	Ø 270 x 235	Ø 13 x 660	DH.WHM121608
for Ki./Py.	4,000	160 - 162	700	170	Ø 325 x 270	Ø 13 x 660	DH.WHM121609
for DU.	5,000	169 - 171	800	170	Ø 325 x 270	Ø 13 x 660	DH.WHM121610

Remarks: DU. = Duran®, Ki. = Kimax, Iw. = Iwaki, Py. = Pyrex

Accessories: Magnetic stirrer bars can be found on page 49