


InLab® Ultra Micro-ISM

pH Ref °C ISM



Micro sensor for very small sample volumes like micro plates, well plates, vials

Scale 1:1



Sensor type combined pH electrode

Order number 30244732

Specifications

pH range 1...11

zero point (pH) 7.00 ± 0.25

temperature range 0...80 °C

temperature probe -

shaft material glass

membrane glass Lot

membrane resistance < 700 M Ω

reference system ARGENTHAL™

reference diaphragm ceramic

reference electrolyte FRISCOLYT-B®

bridge electrolyte -

electrode head MultiPin™

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)

or FRISCOLYT-B® (51340053)