


LE621 IP67

Robust DO sensor for outdoor use

DO °C ISM



Sensor type galvanic DO

Order number 30266975

Specifications

measuring range 0...20 mg/L

zero point (pH) -

temperature range 0...50 °C

temperature probe NTC 30 k Ω

shaft material PVC

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 2.0 m

connector BNC / Cinch

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

10% NaCl