


InLab®605-ISM

dissolved oxygen



Order number 51344611

former order number -

Specifications

measurement range 0 ... 200 %; 0 ... 20 mg/L
compatible with SevenExcellence™ with DO/BOD module
Seven2Go™ S4,
SevenGo Duo pro™
SG68

temperature range 0 ... 60 °C

temperature probe NTC 22 kΩ

shaft material PPS

shaft length 120 mm

shaft diameter 12 mm

connector

electrode head -

cable length 2.0 m

connector BNC/RCA (cinch)

Storage

dry

Special features

IP67

fixed cable

membrane is protected by a steel mesh

shaft made of glass fiber reinforced PPS

compatible with appropriate METTLER TOLEDO meters with Intelligent Sensor Management (ISM) functionality

Choice of accessories

Membrane kit: 51340293 (containing 3 membrane bodies and 25 mL of electrolyte)

Electrolyte: 51340294

Scale 1:1