


# OVERHEAD STIRRERS | Accessories

Product	Ident. No.
<b>TEMPERATURE SENSOR</b>	
H 67.60 stainless steel	0025000215
H 67.61 stainless steel	0020003865
H 68.55 stainless steel	0025001479
H 70 extension cable	0002735600



<b>SENSOR IN COMBINATION WITH H 70</b>	
H 62.51 Stainless steel	0002735451
H 66.51 Glass-coated	0002735551
H 66.53 coated	0004499900



<b>STAND</b>	
<b>Plate stand</b>	
IKA [scale]	0025004318
R 1825	0003160000
R 1826	0003160100
R 1827	0003160200
<b>H-stand</b>	
R 2722	0001412000
<b>Telescopic stand</b>	
R 2723	0001412100
R 474	0001643000
<b>Floor stand</b>	
R 472	0000738700
R 2850	0020002900
R 6547 H	0020018378
R 6547 E 400V	0020018384

MICROSTARS digital   control	MINISTARS digital   control	EUROSTAR 20   40 digital
control	control	-
control	control	-
control	control	-
control	control	-

-	-	-
-	-	-
-	-	-
++	++	++
++	++	++
++	++	++
++	++	++
++	++	++
++	++	++
++	++	++
++	++	++
-	-	-
-	-	-
-	-	-
-	-	-

# OVERHEAD STIRRERS | Accessories

EUROSTAR 60 digital   control	EUROSTAR 100 digital   control	EUROSTAR 200 digital   control	EUROSTAR 400 digital   control	EUROSTAR 100 & 200 control P 4	EUROSTAR 20 high speed digital   control	RW 20 digital	RW 28 digital	RW 47 digital
control	control	control	control	++	control	-	-	-
control	control	control	control	++	control	-	-	-
control	control	control	control	++	control	-	-	-
control	control	control	control	++	control	-	-	-

control	control	control	control	control P 4	20 hs control	-	-	-
control	control	control	control	-	20 hs control	-	-	-
control	control	control	control	control P 4	20 hs control	-	-	-
++	++	++	-	++	++	++	-	-
++	++	++	-	++	++	++	-	-
++	++	++	-	++	++	++	-	-
++	++	-	-	-	++	++	-	-
++	++	++	++	++	++	-	++	-
++	++	++	++	++	++	-	++	-
-	-	-	-	-	-	-	-	++
-	++	++	++	++	-	-	-	-
-	-	-	-	-	-	-	-	++
-	-	-	-	-	-	-	-	++

- Not compatible | + Installable | ++ IKA recommended

<b>CLAMPS // Compatibility with stands (++ / -)</b>					
Product	R 182	R 270	R 271	RH 3	RH 5
IKA scale	++	-	-	++	-
R 1825	++	-	-	++	-
R 1826	++	-	-	++	-
R 1827	++	-	-	++	-
R 2722	-	++	++	-	++
R 2723	-	++	++	-	++

