


Further Accessories



Universal Stand S2

Stand tube Ø 25 mm,
height 700 mm,
leg distance 370 mm,
weight 5.8 kg

P/N 570-12000-00



Stand S2 XXL

Stand tube Ø 25 mm,
height 1.000 mm,
leg distance 370 mm,
weight 6.0 kg

P/N 570-12200-00



Telescope Stand

Stand tube Ø 32 mm,
height 725 mm to 1,025 mm,
leg distance 370 mm,
weight 7.7 kg

P/N 570-12100-00



Clamp

For stand S2, S2 XXL and
telescope stand, Ø 13 – 32 mm

P/N 570-22000-00



Flex Coupling

With clamping spigot, for stirrer
shafts with Ø 10 mm

P/N 509-03000-00



Stirrer Guide (NS 29/32)

For stirrer shafts with Ø 8 mm,
ground PTFE core; suitable for
vacuum, perfect guide for stainless
steel and glass stirrer shafts

P/N 509-09000-00



Shaft Guard

For Hei-TORQUE, made of PMMA,
incl. adapter set, height-adjustable
from approx. 187 – 312 mm

P/N 509-08100-00



Flexible Shaft

Incl. chuck,
1,300 mm overall length

P/N 509-07000-00



RS 232 Cable

9-pin, for Hei-Connect and
Hei-TORQUE Precision models

P/N 14-007-040-72

Adapter set (without illustration)

To fasten the stirrer shaft guard on the Hei-TORQUE Overhead Stirrer

P/N 11-002-501-02