

# Disposable plunging siphon



If strongly adhesive liquids are sampled, it can be a lot less expensive to dispose of the sampler than to carry out a complicated cleaning process.

The disposable dip sampler were developed especially for this purpose. They are made of PE, transparent and the contents are easily visible. PE provides a surface that resists adhesion of the sample material, minimising sample retention on the dip sampler.

## Use

1. Insert the dip sampler into the liquid and fill it.
2. Close the thumbhole above the handle using your thumb.
3. Take out the dip sampler and empty the liquid into the sample bottle.

### Disposable plunging siphon

Volume ml	Immersion depth cm	Total length cm	Tube Ø mm	Quantity scale from	Item no.
75	35	50	27	12	5387-0050 N4
175	60	75	27	12	5387-0075 N4
275	85	100	27	12	5387-0100 N4

