


Sampling bag SteriBag | the sterilised packaging bag




The sterilised sampling bag SteriBag for samples in chemistry, pharmaceuticals, medicine, veterinary science and foodstuffs etc.

The bags are made of polyethylene and are therefore exceptionally tearproof and leakproof.

The SteriBag sampling bag is distinguished by its safe handling. When used according to instructions, samples are guaranteed sterile.

Integrated safety: The tab-ends are sealed in plastic film to prevent skin injuries.

- ▶ Sterile, RNase, DNase and pyrogen-free
- ▶ Can be labelled
- ▶ Leakproof
- ▶ Save space
- ▶ Transparent
- ▶ With practical closure to tear-open
- ▶ Proven millions of times



Sampling bag SteriBag

Capacity ml	Model	L x W mm	Pieces per pack	Item no.	
60	without labelling area	127 x 76	500	5344-0060	N3
60	with labelling area	127 x 76	500	5344-0061	N3
150	without labelling area	178 x 76	500	5344-0180	N3
150	with labelling area	178 x 76	500	5344-0181	N3
450	without labelling area	229 x 114	500	5344-0540	N3
450	with labelling area	229 x 114	500	5344-0541	N3
650	without labelling area	229 x 140	500	5344-0720	N3
650	with labelling area	229 x 140	500	5344-0721	N3
800	without labelling area	305 x 114	500	5344-0810	N3
800	with labelling area	305 x 114	500	5344-0811	N3
1650	without labelling area	305 x 178	250	5344-1650	N3
1650	with labelling area	305 x 178	250	5344-1651	N3