

DIAPHRAGM VACUUM PUMP LABOPORT® N 842.3 FT.18



Technical features

- 100 % oil-free transfer – thereby pure transferring, evacuation and compression of gases
- High level of vapor and condensate compatibility
- PTFE pump head combined with PTFE-coated diaphragm are ideal for extremely aggressive/ corrosive gases and vapors
- Maintenance-free
- Environmentally friendly
- Chemically-resistant
- Gastight, leakage rate approx. 6×10^{-3} mbar x l/s, not tested in serial production

Technical data

| | |
|---|-------|
| Flow rate [m³/h] at atm. pressure | 2.04 |
| Flow rate [l/min] at atm. pressure | 34 |
| Ultimate vacuum [mbar abs.] | 2 |
| Operating pressure [bar g] | 1 |
| Voltage [V] | 230 |
| Frequencies [Hz] | 50 |
| Motor protection | IP 44 |
| Power P ₁ [W] | 245 |
| Operating current [A] | 1.5 |
| With thermal switch and power fuse | |
| Motors with other voltages and frequencies on request | |

MATERIALS

| | |
|-----------|-------------|
| Pump head | PTFE |
| Diaphragm | PTFE-coated |
| Valves | FFPM |

OTHER DATA


| | |
|--|-----------------|
| Weight [kg] | 13.4 |
| Hose connections | ID 10 |
| Dimensions W x H x D [mm] | 167 x 228 x 341 |
| Permissible media and ambient temperature [°C] | +5 ... +40 |

ACCESSORIES

| | ID |
|------------|--------|
| Spares kit | 057359 |


Dimensions

N 842.3 FT.18




Performance characteristics

PUMP DOWN TIME FOR 10 LITER VESSEL



FLOW CAPACITY



YOUR LOCAL PARTNER AROUND THE WORLD

KNF SALES ORGANIZATION

Australia KNF Regional Office
Moreland West VIC 3055
Tel. +61 3 9386 4959
info@knf.com.au
www.knf.com.au

Benelux Netherlands KNF Verder B.V.
3451 GG Vleuten
Tel. +31 30 677 92 40
info@knf-verder.nl
www.knf-verder.nl

Benelux Belgium and Luxemburg
KNF Verder N.V.
2630 Aartselaar
Tel. +32 3 871 96 24
info@knf.be
www.knf.be

China KNF Technology (Shanghai) Co., Ltd.
Shanghai 201203
Tel. +86 21 5109 9695
info@knf.com.cn
www.knf.com.cn

France KNF Neuberger SAS
68128 Village-Neuf
Tel. +33 389 70 35 00
info@knf.fr
www.knf.fr

Germany KNF Neuberger GmbH
79112 Freiburg
Tel. +49 7664 5909 0
info@knf.de
www.knf.de

India KNF Pumps + Systems (India) Pvt. Ltd.
Hinjewadi Pune 411 057
Tel. +91 20 640 13 923
info@knfpumps.in
www.knfpumps.in

Italy KNF Italia S.r.l.
20063 Cernusco s. Naviglio (MI)
Tel. +39 02 929 04 91
info@knf.it
www.knf.it

Japan KNF Japan Co. Ltd.
Tokyo 104-0033
Tel. +81 3 3551 7931
info@knf.co.jp
www.knf.co.jp

Korea KNF Neuberger Ltd.
135-502 Seoul
Tel. +82 2 959 0255
knf@knfkorea.com
www.knfkorea.com

Latin America KNF Regional Office
Tel. +1 609 649 1010
gb@knf.com
www.knf.com/es

Morocco, Tunisia, Algeria
KNF Neuberger SAS
68128 Village-Neuf
Tel. +33 389 70 35 00
info@knf.fr
www.knf.fr

Singapore KNF Regional Office
Tel. +65 9722 1994
info@knf.com.sg
www.knf.com.sg

Sweden, Finland, Denmark, Norway
KNF Neuberger AB
117 43 Stockholm
Tel. +46 8 744 51 13
info@knf.se
www.knf.se

Switzerland KNF Neuberger AG
8362 Balzerswil
Tel. +41 71 973 99 30
knf@knf.ch
www.knf.ch

UK, Ireland KNF Neuberger U.K., Ltd.
Witney, Oxfordshire OX28 4FA
Tel. +44 1993 77 83 73
info@knf.co.uk
www.knf.co.uk

USA, Canada KNF Neuberger, Inc.
Trenton, NJ 08691-1810
Tel. +1 609 890 8600
knfusa@knf.com
www.knfusa.com

KNF PRODUCT CENTERS

Gas Pumps KNF Neuberger GmbH
DE-79112 Freiburg
info@knf.de
www.knf.de

Micro Gas Pumps KNF Micro AG
CH-6260 Reiden
info@knf-micro.ch
www.knf-micro.ch

Liquid Pumps KNF Flodos AG
CH-6210 Sursee
info@knf-flodos.ch
www.knf-flodos.ch