

Spare parts & Accessories

Filling tubes made of FEP

Length x Ø mm	suitable for	Order number
Telescopic 170-330 x 4	NEW LABmax 2.5-10 ml	5 375 000 TK
Telescopic 170-330 x 6	NEW LABmax 25-100 ml, TITREX	5 375 001 TK
310 x 5	LABmax 2.5-10 ml	5 375 000 N 310
390 x 5	LABmax 2.5-10 ml	5 375 000 N 390
310 x 7	LABmax 25-100 ml, TITREX	5 375 001 310
390 x 7	LABmax 25-100 ml, TITREX	5 375 001 390
	MINISPENSOR 50-100 µl	5 375 004 M
	MINISPENSOR 100-2000 µl	5 375 014



Thread adapters made of PP

Thread GL	suitable for	Order number
A22	LABmax, TITREX, MINISPENSOR	5 377 022
A25	LABmax, TITREX, MINISPENSOR	5 377 025
A28	LABmax, TITREX, MINISPENSOR	5 377 028
A30	LABmax, TITREX, MINISPENSOR	5 377 030
A32 / A28*	MINISPENSOR	5 377 032
A38	LABmax, TITREX, MINISPENSOR	5 377 038
A40	LABmax, TITREX, MINISPENSOR	5 377 040
A45 / A32**	LABmax, TITREX	5 497 600
A45	LABmax, TITREX, MINISPENSOR	5 377 045



* from A28 to 32 (MINISPENSOR)
** from A45 to A32 inverted

Additional accessories

Check valve with plastic coating

For fuming and smelling liquids like linoleum, bromine etc.

suitable for	Order number
LABmax, TITREX	5 377 300



Calcium chloride tube with plastic coating

Protection against moisture from the outside

suitable for	Order number
LABmax, TITREX	5 377 310*



One-way filter made of PP,

With membrane made of regenerated cellulose (dust and bacteria filter)

Pore size	Filter Ø	Packaging size pieces	suitable for	Order number
0.2 µm	13 mm	100	LABmax TITREX	8 212 001
0.2 µm	30 mm	100		8 212 002
0.45 µm	13 mm	100		8 212 003
0.45 µm	30 mm	100		8 212 004
0.45 µm	30 mm	50 individual sterile packed		8 214 001
0.2 µm	30 mm	50 individual sterile packed		8 214 002



* delivered without calcium chloride