


# Heating mantle | multi-place

multi-places (3 or 6), individual controller for each place, up to 450°C



## WHM multi-place, 100°C - 450°C

### Ideal for:

- repetitive extraction (Kjeldahl, Soxhlet and other extraction procedures), laboratory reflux and distillation applications
- Duran®, Kimax, Pyrex and Iwaki-brand round flask

### Features:

- nickel chrome heating element
- coil type: acid-resistant, long life time
- highly processed aluminium case offers robust standing
- easy handling
- lower profile for space-saving convenience
- three support rods and two hooks included
- CE certified and unique serial number for tracing

### Controller:

- built-in analog controller



6-places heating mantle  
with three rods and two hooks (included)



3-places heating mantle

### 6-places heating mantle

Flask capacity (ml)	Flask		Watts (W)	Inside depth (mm)	Dimensions W x D x H (mm)	Order number	
	max. Ø (cm)	Bottom				with 230V	with 120V
250/300	83 - 86	flat	840	50	800 x 240 x 110	DH.WHM12291	DH.WHM12292
250/300	83 - 86	round	840	50	800 x 240 x 110	DH.WHM12293	DH.WHM12294
500	102 - 105	round	1440	65	950 x 270 x 130	DH.WHM12295	DH.WHM12296

### 3-places heating mantle

Flask capacity (ml)	Flask		Watts (W)	Inside depth (mm)	Dimensions W x D x H (mm)	Order number	
	max. Ø (cm)	Bottom				with 230V	with 120V
500	102 x 105	round	720	65	550 x 270 x 130	DH.WHM12391	DH.WHM12392
1,000	130 x 131	round	1050	75	690 x 340 x 170	DH.WHM12393	DH.WHM12394

Please note: No clamps, bossheads or extraction apparatus included!

Accessories: Suitable extraction apparatus can be found on page 235

Suitable clamps and bossheads can be found on page 217. Recommended: bosshead 7 400 502 and clamp 7 400 720.