



designed for scientists



ALMEMO 2490

/// Data Sheet

The ALMEMO 2490-2 allows to connect up to two sensors with the special ALMEMO plugs (see list of accessories) for pH and/or Temperature measurements. It displays the readings on the device itself or can be connected to a PC, using labworldsoft. The device is provided with a carrying case, protective holster, and power adapter (90-260 VAC Input / 12 VDC Output 2,5A - Europlug and adapter to US plug).



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	83 x 42 x 127
Weight [kg]	0.26
Permissible ambient temperature [°C]	-10 - +50
Permissible relative humidity [%]	90
RS 232 interface	yes
USB interface	yes

