

## Hei-VAP Ultimate Control



Adds numerous smart features to the Hei-VAP Ultimate, as well as an overview of all parameters and the current process

- 7" touch display for controlling the entire system including the vacuum and the cooling temperature
- For precise monitoring, the pictogram on the display also symbolizes active functions
- Separate knobs for quick access: the right one for direct vacuum control, the left one for the rotation
- Setting of values either via the touch display or the separate quick access knobs
- Direct access to two favorites on the start screen
- Expandable solvent library. You can also transfer your own favorites to other devices using a USB stick
- Define and save own evaporation processes (ramps)
- Automatically evaporate with the Dynamic AUTO<sub>accurate</sub> Program (DAA). Convenient and customizable: from foam avoidance to highly efficient evaporation
- USB and Micro SD plus LAN / RS 232 interfaces (Lab 4.0 ready)

## Hei-VAP Ultimate Control - Technical Data

|  |  |
|--|--|
| Rotation speed indicator                           | 7" LCD color graphic display   |
| Permissible ambient conditions                     | 5 – 31 °C at 80 % rel. humidity<br>32 – 40 °C decreasing lineary up to max 50% rel. humidity |
| Protection class IEC 60529                         | IP 20  |
| Lift   | Hand lift or Motor lift  |
| Height adjustment speed                            | Hand lift models: manual<br>Motor lift models: 25 mm/s                                       |
| Height adjustment                                  | 155 mm   |
| Drive  | brushless DC motor with electronic speed control   |
| Heating capacity                                   | 1,300 W  |
| Temp. range heating bath                           | 20 – 210 °C  |
| Temperature accuracy bath                          | ±1K  |
| Overheat protection bath                           | cut-off at 5 °C over set temperature via separate PT 1000                                    |
| Bath temperature setting                           | 7" LCD color graphic display   |
| Heating bath temperature control                   | electronic / digital   |
| Secondary overtemp. cut off                        | 250 °C   |
| Material heating bath                              | stainless steel V4A  |
| Diameter heating bath                              | 255 mm   |
| Volume heating bath                                | 4.5 l  |
| Display vapor temperature                          | Yes  |
| Required accessories for display vapor temperature | Vapor temperature sensor   |
| Integrated vacuum controller                       | Yes  |
| Timer  | Yes  |
| Programmable ramps                                 | Yes  |
| Supply Power                                       | 1,400 W  |

|   |  |
|---|--|
| Weight without glassware assembly                                   | 15 kg  |
| Dimensions with glassware set G3 (w/d/h)                            | HL: 739 x 532 x 887 mm<br>ML: 739 x 477 x 887 mm |
| Dimensions with glassware set G3 XL (w/d/h)                         | HL: 739 x 532 x 920 mm<br>ML: 739 x 477 x 920 mm |
| Condensing surface  | 1,400 cm <sup>2</sup>                            |
| Condensing surface XL condenser                                     | 2,200 cm <sup>2</sup>                            |
| Max. size evaporating flask   | 5 l  |
| Rotation speed range  | 10 - 280 rpm                                     |
| Integrated Hei-CHILL controller                                     | Yes  |
| Vacuum program Dynamic AUTOaccurate<br>requires AUTOaccurate sensor | Yes  |
| Solvent library   | Yes  |