

Reax control

P/N: 541-11000-00



The precise one with speed scaling

All the outstanding features of the Reax Top's "little brother" can be found in this model: a rotational speed of up to 2,500 rpm and a 5 mm shaking orbit.

The electronic speed control and the scaling of the speed display allow precise settings. Applications that require a determined speed are thus shaken more gently. Shear force-sensitive organisms remain particularly protected by the constant speed in the low range.

Reax control - Technical Data

| | |
|--------------------------------|--|
| Overheat protection | self-resetting |
| Permissible ambient conditions | 5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity |
| Weight | 2.8 kg |
| Protection class IEC 60529 | IP 22 |
| Timer | - |
| Rotation speed range | 0 - 2,500 rpm |
| Power input | 51 W |
| Dimensions (w/d/h) | 134 x 172 x 105 mm |
| Max. Load | - |
| Operating Mode | automatic or continuous |
| Rotation speed setting | analog |
| Orbit / Stroke | 5 mm |
| Platform size | - |
| Accessories included | - |
| Angle | - |
| Motion | vibrating |