

Filtration | glass

Filter Ø25mm, 15ml glass filter holder

Ideal for:


- small volume vacuum filtration for biological or particulate contamination analysis



VF5 - silicone stopper




VF8 - silicone stopper
with stainless steel support screen and PTFE gasket



- 1
- 2
- 3
- 4
- 5

VF5 - 15ml funnel					
#	Item	Material	Size	Order number	
	25mm Glass microanalysis filtration set, complete				00829 167200-05
	25mm Glass microanalysis filtration set, filter holder				00829 167120-32
	Filter holder	Borosilicate glass	Ø25mm		
1	Funnel	Borosilicate glass	15 ml	00829 167220-15	
2	Clamp	Anodized aluminium		00829 167240-05	
3	Support base	Borosilicate glass		00829 167230-04	
4	Stopper	Silicone	No. 5	00829 167110-05	
5	Receiver flask	Borosilicate glass	125 ml	00829 167250-05	

VF8 - 15ml funnel with support screen					
#	Item	Material	Size	Order number	
	25mm Glass microanalysis filtration set, complete				00829 167200-08
	25mm Glass microanalysis filtration set, filter holder				00829 167120-34
	Filter holder	Borosilicate glass	Ø25mm		
1	Funnel	Borosilicate glass	15 ml	00829 167220-15	
2	Clamp	Anodized aluminium		00829 167240-05	
3	Support screen	Stainless steel		00829 167230-32	
4	Gasket	PTFE		00829 167230-42	
5	Support base	Borosilicate glass		00829 167230-08	
6	Stopper	Silicone	No. 5	00829 167110-05	
7	Receiver flask	Borosilicate glass	125 ml	00829 167250-05	



- 1
- 2
- 3
- 4
- 5
- 6
- 7

VF8

Filter Ø47mm, 300ml glass filter holder

Ideal for:


- solvent purification and general filtration for microbiology, suspended solid analysis etc.




VF6 - silicone stopper



VF7 - silicone stopper with stainless steel support screen and PTFE gasket



VF6



VF7

VF6 - 300ml funnel

#	Item	Material	Size	Order number
47mm Glass filtration set, complete				00829 167200-06
47mm Glass filtration set, filter holder				00829 167120-31
1	Filter holder	Borosilicate glass	Ø47mm	
2	Funnel	Borosilicate glass	300 ml	00829 167210-03
3	Clamp	Anodized aluminium		00829 167240-01
4	Support base	Borosilicate glass		00829 167230-03
5	Stopper	Silicone	No. 8	00829 167110-16
6	Receiver flask	Borosilicate glass	1000 ml	00829 167250-03
	Fixing sucker	Silicone		00829 167200-45

VF7 - 300ml funnel with support screen

#	Item	Material	Size	Order number
47mm Glass filtration set, complete				00829 167200-07
47mm Glass filtration set, filter holder				00829 167120-33
1	Filter holder	Borosilicate glass	Ø47mm	
2	Funnel	Borosilicate glass	300 ml	00829 167210-03
3	Clamp	Anodized aluminium		00829 167240-01
4	Support screen	Stainless steel		00829 167230-31
5	Gasket	PTFE		00829 167230-41
6	Support base	Borosilicate glass		00829 167230-07
7	Stopper	Silicone	No. 8	00829 167110-16
8	Receiver flask	Borosilicate glass	1000 ml	00829 167250-03
	Fixing sucker	Silicone		00829 167200-45

Filtration | glass


Filter Ø47mm, 300ml glass filter holder

Ideal for:

- purification of organic, corrosive liquid such as solvent (mobile phase) for HPLC analysis



VF3 - ground joint socket
Ø47mm, 300ml



VF3

VF3 - 300ml funnel					
#	Item	Material	Size	Order number	
	47mm Glass filtration set, complete				00829 167200-03
	47mm Glass filtration set, filter holder				00829 167120-30
	Filter holder	Borosilicate glass	Ø47mm		
1	Funnel	Borosilicate glass	300 ml	00829 167210-03	
2	Clamp	Anodized aluminium		00829 167240-01	
3	Suction cap	Borosilicate glass		00829 167230-01	
4	Receiver flask	Borosilicate glass	1000 ml	00829 167250-01	
5	Fixing sucker	Silicone		00829 167200-45	