

KERN DAB-A01.

KERN

Temperature calibration set consists of measuring sensor and display device



Dimensions packaging (W×D×H)	175×110×55 mm
Net weight	0,215 kg
Shipping method	Parcel service
Net weight approx.	0,25 kg
Gross weight approx.	0,30 kg
Shipping weight	0,26 kg

Pictograms

STANDARD



Category

Brand	KERN
Product category	Balance
Product group	Temperature calibration set

Measuring System

Temperature [Min] (°C)	-200 °C
Temperature [Max] (°C)	1370 °C
Readability temperature [d] (°C)	0,1 °C
Units	°C °F

Approval

CE mark	✓
---------	---

Construction

Dimension (W×D×H)	145×71×35 mm
Material	plastic
Material housing	ABS plastic
Cable length	0,9 m

Power Supply

Supplied power supply	Battery
Battery	4×1.5 V AAA
Battery / accumulator type	Alkali(-Manganese)
Battery connection	Pad + Spring
Battery capacity	1.200 mAh
Battery voltage	6 V
Battery operating time	400 h

Packing & Shipping

Readability force [d] (N)	1 d
---------------------------	-----