

KERN DLB-A01N.

KERN

Temperature calibration set consists of measuring sensor and display device



Category

Brand	KERN
Product category	Balance
Product group	Temperature calibration set

Measuring System

Temperature [Min] (°C)	40 °C
Temperature [Max] (°C)	200 °C
Readability temperature [d] (°C)	1 °C
Units	°C

Construction

Dimension (W×D×H)	100×0×30 mm
Material housing	metal, lacquered
Cable length	0 m

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	215×200×60 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,30 kg
Shipping weight	0,5 kg

Pictograms

STANDARD

