


Product Information



Elmasonic Accessories - Modular basket system

Modular basket

- stainless steel with plastic handles (10 - 120 & 180 - 300) or Rilsan coating (150)
- ergonomic handles
- autoclavable
- Bayonet holes in the bottom for mounting holders (e.g. Elma Pin, silicone holder bottom comb)
- Square openings in the side walls for mounting silicone holders





Basket size	10	20	30	40	60	80	100	120	150	180	300	
Order number	1117923	1119922	1108916	1113045	1113050	1111838	1113387	1113390	1113110	1113404	1113072	
Int. dim. of basket W / D / H (mm)	155 / 70 / 35	110 / 105 / 50	195 / 105 / 50	195 / 105 / 95	255 / 105 / 95	455 / 105 / 55	250 / 195 / 90	250 / 195 / 135	450 / 255 / 55	280 / 255 / 135	450 / 255 / 135	
Chamfer V x W (mm)	--	20 x 20					25 x 25					
Hole diameter bottom Y x Z (mm)	11 x 15											
Hole diameter sides Y ₁ x Z ₁ (mm)	--	10 x 10										
Loading weight max. (kg)	1,0	1,0	1,0	3,0	5,0	5,0	6,0	7,0	7,0	8,0	10,0	
Suitable for unit / bath size	10	15	20	30	40	60	80	100	120	150	180	300
• Elmasonic Easy	•		•	•	•	•		•	•		•	•
• Elmasonic Select				•	•	•		•	•	•	•	•
• Elmasonic P				•		•	•	•		•		•
• Elmadry				•				•				•
• Elma Sink					•	•				•		•

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system

Modular basket

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.