


Product Information



Elmasonic Accessories - Covers

Plastic lid

- available in cobalt blue (RAL 5013)*¹, ultramarine blue (RAL 5002)*² and silver grey (RAL 7001)*³
- reduce the noise of the ultrasonic cleaning
- save energy and costs by collecting and recirculating evaporating water from the ultrasonic bath
- serve upside down as a drip tray for the basket and cleaned parts




	10	15	20	30	40	60	70	80	100	120	150	300	180	200
Order number cobalt blue (RAL 5013)* ¹				1101126		1101129		1110200	1101133		1101135		1101141	1106428
Order number ultramarine blue (RAL 5002)* ²	1003278	1003280		1003281		1003282		1003283	1003285		1003287		1003286	--
Order number silver grey (RAL 7001)* ³	1071693	1071694		1071695		1071696		1110202	1071697		1071699		1071698	--
Ext. dimensions W / D / H (mm)	235 / 115 / 20	220 / 175 / 20		305 / 175 / 20		270 / 190 / 20		570 / 170 / 20	380 / 280 / 20		575 / 335 / 20		405 / 340 / 20	600 / 365 / 20
Suitable for	10	15	20	30	40	60	70	80	100	120	150	300	180	200
• Elmasonic Easy	•		•	•		•			•	•		•	•	
• Elmasonic Select				•	•	•		•	•	•	•	•	•	
• Elmasonic P				•		•	•		•	•		•	•	

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Covers

Plastic lid

Dimension definition



The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.