

Product Information



Elmasonic Accessories - Baskets

Basket

- made of stainless steel and with Rilsan-coated handles



	70	100 medium	120 medium	100 flat	120 flat
Order number	1004247	1071801		1077166	
Basket int. dim. W / D / H (mm)	465 / 105 / 50	255 / 200 / 45		255 / 200 / 25	
Mesh size Y x Z (mm) / Hole diameter Dia (mm)	9 x 1				
Loading weight max. (kg)	5,0	6,0			
Suitable for	70	100	120	100	120
• Elmasonic Easy		•	•	•	•
• Elmasonic Select		•	•	•	•
• Elmasonic P	•		•		•
• Elmadry TD			•		•


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Baskets

Basket



Dimension definition

Basket internal dimensions



Mesh size

Y = Mesh size; Z = Wire diameter (basket sizes 70, 100 medium / flat, 120 medium / flat)



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.