## **Product Information**



### **Elmasonic Accessories - Covers**

#### Stainless steel lid

- reduce the noise of the ultrasonic cleaning
- save energy and costs by reducing waste heat and collecting / recirculating the water vapour



	30	40	60	100	120	150	180	300	500	900
Order number	111	1930	1111931	1111932		1111933	1111934	1111933	1082096	1082091
External dimensions W / D / H (mm)	265 / 1	60 / 55	325 / 175 / 55	325 / 265 / 55		530 / 325 / 55	355 / 325 / 55	530 / 325 / 55	540 / 370 / 50	650 / 550 / 50
Suitable for	30	40	60	100	120	150	180	300	500	900
Elmasonic Easy	•	•	•	•	•		•	•		
Elmasonic Select	•	•	•	•	•	•	•	•	•	•
• Elmasonic P	•		•	•	•		•	•		


# **Product Information**




## **Elmasonic Accessories - Covers**

Stainless steel lid

**Dimension definition** 





The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.