

Product Information



Elmasonic Accessory - Covers



Noise protection covers

- reduce the noise of the ultrasonic cleaning by 40 % (500) and 54% (900)* and retain process vapours
- save energy and costs by reducing waste heat and collecting / recirculating the water vapour
- available as a retrofit kit for self-assembly and space-savings attached to the ultrasonic device



	500	900
Order number	1100349	1100292
External dimensions W / D / H (mm)	580 / 495 / 180	680 / 665 / 180
Height installed H ₁ (mm)	65	
Suitable for	500	900
• Elmasonic Select	•	•

Dimension definition



*Only when closed and without basket / load.

The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.