

Elmasonic Accessories - Modular basket system

Level set

- Level made of stainless steel
- for horizontal division of the modular basket
- Product contents Level set:
 - Level size S and M: 1x Level + 5x Silicone holder per set
 - Level size L and XL: 1x Level + 10x Silicone holder per set
- maximum one horizontal division per modular basket allowed
- the level can be used at different heights
- autoclavable
- Bayonet holes in the bottom for mounting holders (e.g. Elma Pin, silicone holder bottom comb)




Level size	S	M	L	XL
Order number	1125225	1125226	1125227	1125228
Ext. dimension W / D / H (mm)	195 / 105 / 1	255 / 105 / 1	250 / 195 / 1	285 / 250 / 1
Hole diameter Dia (mm)	11 x 15			
Weight per level (g)	110	150	280	410
Max. loading weight per basket size and level (kg)	1,0	1,0	2,0	2,0

Level size	S	M	L	XL
Suitable for basket size:				
30	1 set			
40	1 set			
60		1 set		
80	2 sets side by side			
100			1 set	
120			1 set	
150			2 sets side by side	
180				1 set
300			2 sets side by side	

Elmasonic Accessories - Modular basket system


Level set

Dimension definition




Example of application

Basket size 60 with Level set size M



Basket size 300 with 2 Level sets size L side by side



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.