


Elmasonic Accessories - Modular basket system

Dividing wall set

- Dividing wall made of stainless steel
- for vertical division of the modular basket
- Product contents Dividing wall set: 2x Dividing wall + 10x Silicone holder per set
- autoclavable
- Square openings in the side walls for mounting silicone holders




Dividing wall size	S low	S high	M	L low	L high
Order number	1125220	1125221	1125223	1125222	1125224
Ext. dimension W / D / H (mm)	105 / 50 / 1	105 / 90 / 1	200 / 60 / 1	200 / 90 / 1	255 / 130 / 1
Hole diameter Dia (mm)	10 x 10				
Weight per dividing wall (g)	70	40	130	100	230

Dividing wall size	S low	S high	M	L low	L high
Suitable for basket size:					
20	1 Set				
30	2 Sets				
40	2 Sets	2 Sets			
60	2 Sets	2 Sets			
80	4 Sets	4 Sets			
100	2 Sets	2 Sets	2 Sets		
120	2 Sets	2 Sets	2 Sets		
150	4 Sets			4 Sets	
180	2 Sets	2 Sets			2 Sets
300	4 Sets	4 Sets			2 Sets

Elmasonic Accessories - Modular basket system


Dividing wall set

Dimension definition




Example of application

Basket size 300 with 2x Dividing wall set L high + 1x Dividing wall set S high



Basket size 30 with 2x Dividing wall set S low



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.