


Elmasonic Accessories - Modular basket system

Holder set


- made of stainless steel
- Holder hooks Rilsan-coated
- Holder hooks with 10 hooks
- Holder hooks for holding e.g. jewelry and watch parts
- Holder shaft with 4 shaft
- Holder shaft for holding e.g. Impression tray
- Product contents Holder set: 1x Holder shaft, 1x Holder hooks, 5x Silicone holder per set




	Holder set
Order number	1125219
Ext. dimension Holder shaft W / W / H (mm)	105 / 3 / 15
Ext. dimension Holder hooks W ₁ / D ₁ / H ₁ (mm)	110 / 55 / 25
Weight bracket shaft (g)	7
Weight holder hooks (g)	9

dimension definition

Holder shaft



Holder hooks



The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Modular basket system

Holder set

Example of application

Basket size 60 with 1x Holder set



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.