


Product Information



Elmasonic Accessories - Modular basket system

Elma Pin Set

- made of plastic
- Color: Black
- For mounting and holding various components
- Product contents Elma Pin Set: 10x Elma Pin per set
- Bayonet lock for easy mounting on the bottom of the baskets and on the levels



Elma Pin	0509	0158	0767	0613 + 0614		0741	0772	0021R	0001R
Order number	1125209	1125210	1125211	1125212		1125214	1125213	1125216	1125215
Ext. dimension Dia / H (mm)	30 / 35	50 / 40	85 / 65	55 / 45	75 / 45	25 / 30	25 / 55	25 / 55	40 / 75
int. dimension Dia ₁ / H ₁ (mm)	20 / 20	40 / 25	50 / 35	0,35 / 25	25 / 25	4 / 15	5 / 40	15 / 40	30 / 65
Outer overlay Dia ₂ (mm)	--	--	--	15	40	--	--	--	--
Weight per Elma Pin (g)	3	6	7	6		2	3	4	8
									


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Modular basket system


Elma Pin Set

Dimension definition


Elma Pin 0509 & 0158




Elma Pin 0767




Elma Pin 0613 & 0614



Elma Pin 0741 & 0772



Elma Pin 0021R & 0001R



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.