




Elmasonic Accessories - Modular basket system



Silicone holder set (Bottom comb / Side comb)

- for holding various instruments, containers (e.g. flasks) and workpieces
- Scope of delivery silicone holder set bottom comb: 10x silicone holder bottom comb each per set
- Scope of delivery silicone holder set side comb: 5x silicone holder side comb each per set
- Color: White



dimension definition

	Bottom comb	Side comb
Order number	1125217	1125218
Ext. dimension W / D / H (mm)	105 / 10 / 35	150 / 10 / 45
Weight per silicone holder (g)	35	71
		



The drawings are meant as a visualization aid and are not true to scale.


Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.

Elmasonic Accessories - Modular basket system


Silicone holder set (Bottom comb / Side comb)

Example of application

Basket size 150 with 2x Bottom comb



Basket size 60 with 1x Side comb



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.