

Product Information



Elmasonic Accessories - Holders and product carriers


Equipment holder

- made of stainless steel and with Rilsan-coated handles
- for holding laboratory vessels e.g. glass flasks up to 1000ml



Order number	1049704
Basket int. dimension W / D / H (mm)	195 / 95 / 75
Suitable for	
Elmasonic Select, Easy, P	
60	1 x
120	2 x
180	2 x
300	4 x

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma Units.