

Elmasonic Accessories - Holders and Product Carriers

Hooks holder

- with 10 hooks that can be flexibly attached
- for cleaning e.g. rings, chains, etc.
- plastic support surfaces to prevent scratching of the tray and the goods




Unit size	60	100	120	150	180	300
Order number	1117432					
Ext. dimension W / D / H (mm)	355 / 60 / 75					
Support height H ₁ (mm)	55					
Dimensions hooks Dia / D ₁ (mm)	1,7 / 55					
Max. length with 10 hooks W ₁ (mm)	250					
Suitable for	60	100	120	150	180	300
• Elmasonic Easy	1 x	2 x	2 x	5 x	4 x	5 x
• Elmasonic Select	1 x	2 x	2 x	5 x	4 x	5 x
• Elmasonic P	1 x		2 x		4 x	5 x

Elmasonic Accessories - Holders and Product Carriers

Hooks holder

Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.