

# Product Information



## Elmasonic Accessories - Holders and product carriers

### Test tube holder

- made of stainless steel
- for test tubes  $\varnothing$  16 mm /  $\varnothing$  25 mm




	30	40	60	80	100	120	150	180	200	300
Order number	1111973				1112084		1112123			
External dimensions W / D / H (mm)	200 / 90 / 95				300 / 90 / 95		360 / 90 / 95			
Height hangers H <sub>1</sub> (mm)	70				70					
Diameter hole Dia (mm)	17				17 / 26					
Max. number of test tubes	16 x $\varnothing$ 16 mm				8 x $\varnothing$ 16 mm + 7 x $\varnothing$ 25 mm					
Suitable for	30	40	60	80	100	120	150	180	200	300
• Elmasonic Easy	1 x	1 x	2 x	4 x	2 x	2 x		2 x		4 x
• Elmasonic Select	1 x	1 x	2 x	4 x	2 x	2 x	4 x	2 x		4 x
• Elmasonic P	1 x		2 x	4 x		2 x		2 x		4 x

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.

## Elmasonic Accessories - Holders and product carriers

### Test tube holder

#### Dimension definition



#### Example of application



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.