Product Information



Elmasonic Accessories - Inserts

Hostal tub

- for cleaning of smallest parts
- protects sensitive parts from damage



	Small		Big	Big with long temples
Order number	2070390000		207040000	2070600000
Ext. dimension W / D / H (mm)	185 / 105 / 95		280 / 105 / 100	280 / 105 / 125
Int. dimension W ₁ / D ₁ / H ₁ (mm)	90 / 90 / 55		195 / 90 / 55	
Hole diameter Dia (mm)	2			
Suitable for	15	20	30	40
Elmasonic Easy		•	•	•
Elmasonic Select			•	•
• Elmasonic P			•	

Product Information



Elmasonic Accessories - Inserts

Hostal tub

Dimension definition





