


## Elmasonic Accessories - Practicals


### Cooling coil

- made of stainless steel
- for connection to a closed-circuit recirculating chiller (cryostat)
- for reducing or stabilising the process temperature in the ultrasonic bath
- use of several cooling coils in combination is possible
- optimal for thermostable processes



	40	70	150	60	80	100	120	130	180	300
Order number	1038735			1037515						
External dimensions W / D / H (mm)	180 / 25 / 130			200 / 25 / 175						
Outer diameter Dia (mm)	10									
Suitable for	40	70	150	60	80	100	120	130	180	300
• Elmasonic Easy	•			•		•	•		•	•
• Elmasonic Select	•		•	•	•	•	•		•	•
• Elmasonic P		•		•			•		•	•

### Dimension definition



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.