

Elmasonic Accessories - Inserts

Immersion baskets


- made of plastic and stainless steel
- can - depending on the device size:
 - be hung individually or several times in the ultrasonic bath,
 - be placed individually or several times in a cleaning basket or
 - be immersed individually in a cleaning glass



Order number	2070380000	1117636	2070540000	1116797	1007545
Ext. dim. W / H (mm)	90 / 180	55 / 75	75 / 140	85 / 105	90 / 160
Int. dim. Dia / H ₁ (mm)	70 / 40	45 / 30	55 / 60	65 / 60	70 / 60
Material	plastic	stainless steel			
Diameter Dia ₁ / Mesh size (mm)	2,4	1x1x0,5			
Suitable for cleaning glass	1000 ml (1047376)	150 ml (1111834)	600 ml slim / high (1117632)	600 ml wide / low (1047272)	1000 ml (1047376)

Elmasonic Accessories - Inserts

Dimension definition



Example of application



The drawings are meant as a visualization aid and are not true to scale.

The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.