

## Elmasonic Accessories - Inserts


### BD-Steribox

- made of plastic
- with opening in the lid
- with removable sieve for rinsing the contents
- for use in the ultrasonic bath placed directly in the basket or hung over an insert cover



Order number	1105991
Ext. dimensions Dia / H (mm)	95 / 85
Int. dimensions Dia <sub>1</sub> / H <sub>1</sub> (mm)	75 / 55
Opening dimensions Cover Dia <sub>3</sub> (mm)	10
Int. dimensions sieve Dia <sub>2</sub> / H <sub>2</sub> (mm)	70 / 45
Opening dimensions sieve bottom W / D (mm)	1,8 / 1,8
Capacity max. (ml)	300

### Dimension definition



The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.