

## Elmasonic Accessories - Inserts

### Cleaning glasses



- can be used for simultaneous cleaning of parts with different chemicals, for removing abrasive substances or for cleaning with aggressive chemicals (e.g. acids or alkalis)
- can either be placed directly in a cleaning basket or in combination with an insert cover be hung in the ultrasonic bath
- Plastic lid and rubber ring included



	150 ml (without plastic cover)	600 ml	600 ml	1000 ml
Order number	1111834	1117632	1047272	1047376
Glas with cover dimensions ext. Dia / H / D (mm)	80 / 85 / 95	100 / 155 / 110	110 / 135 / 125	110 / 185 / 125
Glas dimensions internal Dia <sub>1</sub> / H <sub>1</sub> (mm)	55 / 60	75 / 145	80 / 105	85 / 160
Rubber ring external diameter Dia <sub>R</sub> (mm)	70		100	

### Dimension definition

Cleaning glass  
with plastic  
cover



Rubber ring



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.


# Product Information



## Elmasonic Accessories - Inserts

Cleaning glasses

Example of application



The parts shown are for illustration purposes only and are not included in the scope of delivery.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.