


## Elmasonic Accessories - Inserts

### Insert cover


- made of stainless steel
- with 4 rubber rings
- for hanging in the ultrasonic bath
- for 4 glasses or 4 cleaning cups




	100	120	180
Order number	1111909		1111921
External dimensions W / D / H (mm)	350 / 230 / 40		380 / 20 / 35
Diameter cover hole Dia (mm)	95		
External Diameter rubber ring Dia <sub>1</sub> (mm)	100		
Suitable for	100	120	180
• Elmasonic Easy	•	•	•
• Elmasonic Select	•	•	•
• Elmasonic P	•	•	•

### Dimension definition

Insert cover



Rubber ring



The drawings are meant as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma accessories.