

## **LABSOLUTE® SENSITIVE, NITRILE, POWDER FREE**

Article numbers: 7.696 904 (XS), 7.696 900 (S), 7.696 901 (M), 7.696 902 (L),  
7.696 903 (XL)

### **General information:**

Type: disposable gloves / examination gloves  
Description: All information is printed on the dispenser box.  
Shape: ambidextrous, straight fingers  
Material: Nitrile Butadiene Rubber (NBR)  
Color: lavender blue (valid for HACCP Standard)  
Inner: online chlorinated, powder free  
Outer: no treatment  
Cuff/Surface: rolled cuff, textured fingers  
Length: Ø 240 mm  
Thickness (double): Finger min. 0.14 mm, Palm min. 0.12 mm, Cuff min. 0.10 mm  
Gloves width: median (±10 mm): XS <80 mm, S 80 mm, M 95 mm, L 110 mm, XL >110 mm  
Force at break (EN 455-2): Ø 6,0 N (during shelf life)  
Residual powder: ≤ 2 mg / gloves  
Available sizes: XS (5 - 6), S (6 - 7), M (7 - 8), L (8 - 9), XL (9 - 10)

### **Storage:**

Storage life: 3 years (expiry date is printed on the dispenser box).  
Storage conditions: cool (max. 30 °C), dry, dust free, avoid UV radiation, contact with oxidizing agents, do not store close to photocopy equipment.  
Miscellaneous: Copper ions can discolour the gloves.

### **Directives and Certification:**

Directives: Regulation (EU) 2016/425, CE category III  
Standards: EN 420:2003+A1:2009, ISO 374-1:2016+A1:2018 (type B), EN 374-2:2014 (Level 2 = AQL 1,5), EN 16523-1:2015+A1:2018, EN 374-4:2013, EN ISO 374-5:2016, EN455-1, -2, -3  
Food regulation: Regulation (EU) No. 1935/2004

### **Warning:**

Gloves contain the following allergenic substances: Dithiocarbamate, Zinc-Mercaptobenzothiazol  
This product contains allergenic substances which can cause allergic reactions, including anaphylactic responses.